

Scanning & Analysis Order Form

Name	
Institution	
Address	
Country	
Phone Number	
E-Mail address	
Date sent	
Number of Plates sent	
Send invoice to	
VAT-Number	

 Full Scanning & Analysis

 Scanning only

No	Name to be filed	Species	Cytokine	Color system	Position of negative controls
1					
2					
3					
4					
5					

...continue on page 2 for more plates...

Plate numbers which need to be analyzed with identical parameters:

(Necessary if the plates belong to one and the same experiment/series of experiments and need to be compared directly)

Do some plates need to be analyzed with the same parameters as plates from previous shipments?

If so, please specify the **registration number** and file name of the previously analyzed plate(s):

<input type="checkbox"/>	Shipping via regular mail
<input type="checkbox"/>	Shipping via Express Courier
<input type="checkbox"/>	Data can be burned on CD/DVD
<input type="checkbox"/>	Provide data as Web-Download

Date, Signature

Plate Shipment Requirements:

! The plates have to be completely dry, and packages must include the lid and the bottom of the plate. Clearly label the plates with the above assigned plate number. While labeling the plate, make sure not to write too close to the well edge. Do not put sticky tape on the membranes. Send the plate(s) protected from light (in ImmunoProtect Plate Pouches Part # IPST10 or aluminium foil) and wrap them to prevent them from breaking. Please be aware that CTL cannot guarantee the quality of the scanning and analyzing process if these instructions are not followed.

We recommend sending the plate(s) not later than 1-2 months after development, to avoid loss of information due to the fading of spots.

CTL Europe GmbH, Hans-Boeckler-Str. 19, 53225 Bonn, Germany

Tel +49-(0)228 85 44 78-0 · Fax +49-(0)228 85 44 78-22 · www.immunospot.com · ctl-europe@immunospot.com

Name	
Institution	

No	Name to be filed	Species	Cytokine	Color system	Position of negative controls
6					
7					
8					
9					
10					
11					
12					
13					
14					
15					
16					
17					
18					
19					
20					

Plate numbers which need to be analyzed with identical parameters:

(Necessary if the plates belong to one and the same experiment/series of experiments and need to be compared directly)

**Do some plates need to be analyzed with the same parameters as plates from previous shipments?
If so, please specify registration number and file name of the previously analyzed plate(s):**