



Cellular Technology Ltd.

PRODUCT DATA SHEET – CTL-Wash Supplement Medium

| | |
|---------------------------|--|
| Catalog #: | CTLW-010 |
| Product: | CTL-Wash Medium (available in 100ml) |
| Formulation: | CTL Wash™ contains insulin, transferrin, and glucose among other defined proteins, saturated and unsaturated fatty acids and sugars in a proprietary composition to provide a balanced environment for PBMC. |
| Storage: | 4°C (Refrigerated) Protected from light. |
| Recommended usage: | CTL Wash™ is a 10x serum free formulation to be added to PBS. For best results, it is diluted 1 + 9 in warm (20-37°C) PBS and used warm for washing. Diluted CTL Wash™ can be safely stored at 4°C. CTL Wash™ is sold in convenient sterile 100 ml units ready for dilution in 900 ml sterile PBS. |
| Working solution: | Detailed protocols are available from CTL for all steps of work with PBMC under serum free conditions. |
| Warning: | For laboratory research use only, not for diagnostic use! Not for use on humans. |

Cellular Technology Ltd.
20521 Chagrin Blvd.
Shaker Heights, OH 44122
USA